


HighDenseMEA

HD30/10iR-ITO

High Density Microelectrode Array with Internal Reference Electrode



Technical Specifications: HighDenseMEA

Temperature compatibility	0 - 125 °C
Dimension (W x D x H)	49 mm x 49 mm x 1 mm
Thickness (region of electrodes)	50 µm
Base material	Glass
Contact pads and track material	Indium tin oxide (ITO)
Electrode diameter	10 µm
Interelectrode distance (centre to centre)	30 µm
Distance between electrode fields	500 µm
Electrode height	Planar
Electrode type	Titanium nitride (TiN)
Isolation type	Silicone nitride 500 nm (PEVCD)
Electrode impedance	Approximately 250 - 400 kΩ
Electrode layout grid	2 x (5 x 6)
Number of electrodes	60
Reference electrodes	with internal reference electrode

MC_Rack:

Source Layout in Data Source Setup
MCS Channel map


2 dim. (MEA)
HighDenseMEA.cmp
HighDenseMEA_L.cmp
HighDenseMEA_R.cmp

MEA perfusion chamber

(w/o) Without ring
(gr) Glass ring: ID +/- 20 mm, OD 24 mm, height 6 / 12 mm
(pr) Plastic ring without thread: ID 26 mm, OD 30 mm, height 6 / 3 mm
(pr-T) Plastic ring with thread: ID 25 mm, OD 30 mm, height 6 / 15 mm

Electrode Layout

MEA1060 pins **33 21 32 31 44 43 41 42 52 51 53 54 61 62 71 63**
 Electrode # A3L C3L A4L A5L B4L B5L C4L C5L C1R C2R B1R B2R A1R A2R C3R A3R



The first letter of the electrode number code refers to the row number, the digit is the column number, and the second letter refers to the electrode field (left or right) of the HighDenseMEA. The specified MEA1060 pin numbers are the channel numbers that are used in the MC_Rack program. The electrode D1 of the left electrode field, connected to channel 15 in MC_Rack is missing. It is replaced by a big internal reference electrode.