



HexaMEA40/10

HexaMEA40/10iR-ITO

Hexa Microelectrode Array



Technical Specifications HexaMEA40/10

Temperature compatibility

0 - 125 °C

Dimension (W x D x H)

49 mm x 49 mm x 1 mm

Base material

Glass

Contact pads and track material

ITO (Indium tin oxide)

Electrode diameter

10 µm

Interelectrode distance (centre to centre)

40 µm

Electrode height

Planar

Electrode type

TiN (Titanium nitride)

Isolation type

Silicon nitride 500 nm (PEVCD)

Electrode impedance

Approximately 250 - 400 kΩ

Electrode layout grid

Hexagonal

Number of electrodes

60

Reference electrode

with internal reference electrode (iR)

MC_Rack:

Source layout in "Data Source Setup"

2 dim (MEA)

Channel map

HexaMEA40/10.cmp

MEA perfusion chamber

(w/o) Without ring

(gr) Glass ring: ID +/- 20 mm, OD 24 mm, height 6 / 12 mm

(pr) Plastic ring without thread: ID 26 mm, OD 30 mm, height 6 / 3 mm


(pr-T) Plastic ring with thread: ID 25 mm, OD 30 mm, height 6 / 15 mm

HexaMEA40/10

HexaMEA40/10iR-ITO

Electrode Layout

MEA1060 pins (33) (21) (32) (31) (44) (43) (41) (42) (52) (51) (53) (54) (61) (62) (71) (63)
Electrode # C3 B2 A1 D5 A2 B3 C4 A3 B4 A4 C5 B5 A5 E4 B6 C6



The letter-digit code is the electrode identifier and refers to the position of the electrode in the hexa grid.
The specified MEA1060 amplifier pin numbers are the channel numbers that are used in MC_Rack.